



**Technical Note 179** 

## Fluorescein Performance Data

### Introduction

Fluorescein and its derivative fluorescein isothiocyanate (FITC) are ubiquitous in cell and molecular biology as a fluorescent label for cells and in many other applications. Fluorescein is a fluorescent compound that is easily detected using a fluorometer or fluorescence microscope. Fluorescein has an excitation maximum at 490 nm and an emission maximum at 514 nm.

DeNovix QFX Fluorometer and DS-11 Series Spectrophotometer / Fluorometers enable precise fluorescein quantification using a proprietary optical fluorescence core and four user-selectable LED channels. The data presented in this note was measured in the Basic Fluorometer app of EasyApps<sup>™</sup> software.

### **Materials and Methods**

A fluorescein dilution series diluted from a fluorescein standard (ThermoFisher Scientific cat #F36915) was quantified on a DeNovix QFX fluorometer.

The 50  $\mu$ M standard was diluted in 0.05 M borate buffer with pH 8.5. A dilution series with 12 fluorescein concentrations between 0 nM and 2500 nM was created by serial dilution. Sample volumes of 200  $\mu$ L were measured in 0.5 mL PCR tubes (DeNovix cat #TUBE-PCR-0.5-500). Measurements were made in the Basic Fluorometer app using the blue channel. The blue channel excitation range is 442 – 497 nm, and the emission range is 514 – 561 nm.

### Results

The fluorescein standard curve from 0 to 500 nM is shown in Figure 1. Linearity was exhibited between 0.05 and 2500 nM fluorescein. Table 1 shows the raw data for the full dilution series.

Table 1: Fluorescein Raw Data through Full Concentration Range

| Fluorescein Concentration (nM) | RFU (n = 3) | StDev   | %CV |
|--------------------------------|-------------|---------|-----|
| 0                              | 116.67      | 0.80    | 0.7 |
| 0.05                           | 123.29      | 0.99    | 8.0 |
| 0.1                            | 131.21      | 0.38    | 0.3 |
| 0.25                           | 155.95      | 0.62    | 0.4 |
| 0.5                            | 196.84      | 0.33    | 0.2 |
| 1.0                            | 274.67      | 1.23    | 0.4 |
| 5                              | 901.25      | 1.90    | 0.2 |
| 10                             | 1699.35     | 3.65    | 0.2 |
| 25                             | 4087.39     | 15.74   | 0.4 |
| 100                            | 23394.95    | 33.14   | 0.1 |
| 500                            | 113476.45   | 190.03  | 0.2 |
| 1000                           | 221854.98   | 507.98  | 0.2 |
| 2500                           | 534489.98   | 1459.92 | 0.3 |

# Fluorescein Linearity Full Concentration Range




Figure 1: Fluorescein linearity through the concentration range of 0-500 nM.

29-APR-2022