


Topics for Realizing Robust & Efficient Cell Manufacturing Processes


- Principles of cell therapy manufacturing process optimization/scale-up and a roadmap toward developing robust/efficient processes
- Scale-up challenges for large scale processes
 - Equipment, Materials and Protocols** – Evaluating and selecting high-performing, truly scalable technologies (e.g., bioreactor) and establishing robust and efficient protocols (including seed train culture processes)
 - Hydrodynamics** – Maintaining homogeneous hydrodynamic environment with uniform distribution of turbulent energy dissipation and consistently minimal shear forces across various scales
 - Gas Supply** – Understanding culture scale and achievable cell density via headspace gassing, and finding a new approach for proper gas supply (oxygen) / stripping (carbon dioxide) at large scale in which headspace gassing is not sufficient to maintain oxygen and CO₂ at desired levels
 - Optimal Protocols of Medium Exchange and Cell Harvest** – Finding a new way for rapid and efficient medium exchange (MX) and harvest for cells that are sensitive to process time – e.g., induced pluripotent stem cells (iPSCs) growing as aggregates that tend to fuse with others to result in larger aggregates when they are settled on the bottom of a bioreactor during MX or harvest.

Process Optimization to Develop a Robust & Efficient Process


Importance of Process Development Strategies: Example Case Studies with Different Outcomes



Principle of Process Development




Process Development Roadmap



A Case Study: hMSC Process Optimization and Scale-Up

Figure 4. A summary of a case study of process optimization and scale-up carried out and PIVs investigated to develop a clinical-scale (>10 B cells/lot) manufacturing of hMSCs, in comparison with a process designed based on data from a separate study.



Scale-Up Challenges and Potential Solutions

Selecting a Scalable Bioreactor System (with Similar Geometries) and Determining Agitation Speed to Offer Comparable Hydrodynamics





Figure 6. Observed morphologies of iPSC aggregates for different combinations of volume and agitation rate. Uniformly spherical aggregates correspond to VA EDRs that fall within the optimal range. There is also an inverse correlation between average EDR and aggregate size. Photomicrographs were taken at 10X magnification. Scale bars: 100 μ m.

Strategies for Optimal Gas Supply for Target Culture Scales and Cell Densities

Understanding culture scales & achievable cell densities via headspace gassing: oxygen transfer rate (OTR) vs uptake rate (OUR)



Experimental Condition:

- Vessel: PBS-3 MAG bioreactor
- Working Volume: 3L
- Medium: Water
- Gas inlet: Pure oxygen via headspace (0.5 L/min)
- Agitation speed: 18-45 RPM


Assumptions:

- Specific oxygen consumption rate, q_{O2}: 2E-10 mmol O₂/cell/hr
- Oxygen concentration at 100% DO: 0.2 mmol O₂/L
- Max. oxygen concentration with pure oxygen via headspace: 1.0 mmol O₂/L

Figure 7. kLa values measured in a 3L VW bioreactor under different agitation speeds.

Table 1. Theoretical cell densities achievable under the given experimental conditions and assumptions.

Oxygenation of culture media outside a large bioreactor using an oxygenator



Assumptions:


- Specific oxygen consumption rate, q_{O2}: 2E-10 mmol O₂/cell/hr
- Oxygen concentration at 50% DO set point, C_{O2}: 0.1 mmol O₂/L
- Oxygen concentration in medium oxygenated (w/ pure oxygen in the oxygenator) entering into bioreactor, C_{O2, in}: 1.0 mmol O₂/L

X _i cell density (cells/mL)	F/V (hr ⁻¹)	Flow rate required for different bioreactor working volume (mL/sec)			
		3L	15 L	80L	
500,000	0.11	0.09	0.46	2.47	
750,000	0.17	0.14	0.69	3.70	
1,000,000	0.22	0.19	0.93	4.94	
1,250,000	0.28	0.23	1.16	6.17	
1,500,000	0.33	0.28	1.39	7.41	
1,750,000	0.39	0.32	1.62	8.64	
2,000,000	0.44	0.37	1.85	9.88	
2,250,000	0.50	0.42	2.08	11.11	
2,500,000	0.56	0.46	2.31	12.35	

Table 2. Medium flow rate calculated for different bioreactor working volume for its sufficient oxygenation

Developing Optimal Protocols for Medium Exchange and Cell Harvest for Very Sensitive Cells

Current protocol using Adjustable Height Harvest Valve (AHHV) for rapid and efficient medium exchange and cell harvest in single-use VW bioreactors



Rapid separation of cells from medium outside a bioreactor to enable rapid and complete medium exchange and return to a 2nd bioreactor that contained pre-conditioned fresh medium




Figure 10. A methodology for using an external separation and retention device in conjunction with multiple bioreactors to facilitate rapid, scalable medium exchange

ACKNOWLEDGMENTS

CONCLUSIONS

In order to reliably provide PSC-derived allogeneic cell therapies to vast numbers of patients, a series of optimized unit operations at various scales will be needed to meet target manufacturing lot sizes. Numerous manufacturing processes such as cell aggregate expansion and differentiation, medium and gas exchanges, and cell harvesting all need to be developed and optimized for large-scale use. The proper combination of single-use bioreactor technology and methodologies can avoid potential upstream process bottlenecks and enable robust commercial-scale manufacturing of therapeutic cells.