


## Using CFD Modeling as a Scale-up Tool


The therapeutic use of stem cells has become increasingly popular as a possible treatment option for a variety of diseases and injuries. Typically,  $10^6$  to  $10^9$  cells are required per patient, per treatment [1]. To achieve these numbers, large-scale production techniques must be implemented. Cellular aggregates can be expanded in stirred suspension bioreactors, achieving significantly greater cell quantities per volume compared to traditional static culture [2]. A change in bioreactor scale, however, results in altered hydrodynamic properties that affect cell products. Shear from the local fluid velocity gradient and cell interactions with turbulent eddies influence aggregate sizes, which impact cell pluripotency, differentiation, and proliferation potential [3]. Current bioreactor scale-up equations used for predicting agitation rates maintain single quantity variables that do not capture the variation and complexity of the bioreactor environment. As a result, predicted scale-up agitation rates are often ineffective and damage the cells. Computational Fluid Dynamic (CFD) modeling allows the user to customize geometry so that scale-up equations can be derived between reactors of any given shape and size. We have recently published data [4] that suggests maintaining the volume average energy dissipation rate, derived from CFD simulations, provides a robust method for scale-up of aggregate culture.



## Project Overview


In this study:

- Ansys Fluent software was used to develop a complete (20rpm-100rpm) CFD model of the 0.1L PBS vertical-wheel reactor.
- Model variables (volume average and distributions) were compared to results from the 100mL NDS and the 500mL DasGip bioreactors.
- A hiPSC aggregate pre-formation protocol was developed (comparing pre-formation time and cell density) for bioreactor culture.
- hiPSCs were cultured in 0.1L PBS vertical-wheel reactors at various agitation rates to compare growth and aggregate size.




## Single Use Vertical-Wheel Reactor Model (0.1L)


### Methods



### Velocity Profile




### Shear Stress and Energy Dissipation Rate




## Bioreactor Model Comparison



### Bioreactors Modeled



In this study, the 0.1L PBS vertical-wheel bioreactor was modeled at various agitation rates using CFD software Fluent. Results were compared to models generated for traditional horizontal-blade impeller 100mL NDS and 500mL DasGip bioreactors (modeled with the addition of control probes). It was discovered that while the volume average shear rates between the modeled reactors remained similar, the volume average energy dissipation rate in the 500mL DasGip was nearly 10 times higher than in the 100mL reactors. This increase can likely be attributed to the addition of control probes that occupy a large portion of the liquid headspace in the DasGip. These results suggest that a drastic decrease in agitation rate between the 100mL and 500mL reactors would be required for aggregate size to be maintained throughout scale-up.



## Equation Comparison and Hydrodynamic Variable Distributions




## Discussion


Data generated from the CFD model of the 0.1L PBS Biotech vertical-wheel reactor was used to gain a detailed understanding of the hydrodynamics within the reactor. Unlike traditional horizontal-blade impellers, streamlines in the vertical-wheel reactor form a figure-8 pattern, moving through the entire reactor height with each rotation. The volume average energy dissipation rate (a key variable controlling aggregate size) between the horizontal-blade and vertical-wheel 100mL reactors where found to be very similar. In contrast, the 500mL reactor modeled with control probes, had a much greater volume average energy dissipation rate and should be operated at an agitation rate of approximately half that of the 100mL reactors to maintain aggregate size throughout scale-up. Utilizing our understanding of the bioreactor hydrodynamic forces and control of environmental factors (nutrient and oxygen availability) we will continue to optimize growth and scale-up conditions for hiPSCs cultured as aggregates in the PBS Biotech vertical-wheel bioreactors at various scales.

## Growth of hiPSCs

### Pre-forming Aggregates – Time and Density Results



## Affects of Hydrodynamics on Growth and Aggregate Size



## References

- [1] E.S. Kehoe, D.E., Jing, D., Lock, L.T., Tzanakakis, Scalable stirred-suspension bioreactor culture of human pluripotent stem cells, Tissue Engineering Part A. 16 (2010) 405.
- [2] C.A.V. Rodrigues, T.G. Fernandes, M.M. Diogo, C.L. da Silva, J.M.S. Cabral, Stem cell cultivation in bioreactors, Biotechnology Advances 29 (2011) 815–829.
- [3] S. Stolberg, K.E. McCloskey, Can shear stress direct stem cell fate, Biotechnology Progress. 25 (2009) 10–9.
- [4] B.S. Borys, E.L. Roberts, A. Le, M.S. Kallos, Scale-up of embryonic stem cell aggregate stirred suspension bioreactor culture enabled by computational fluid dynamics modeling, Biochemical Engineering Journal. 133 (2018) 157-167.

## Acknowledgments

