



## Technical Note 181

## DeNovix Trypan Blue Assay Protocol

### Introduction

The Trypan Blue exclusion assay distinguishes between live (unstained) and dead (stained) cells and enables a viability assessment of a cell suspension. Trypan Blue permeates the compromised membranes of dead cells and binds to intracellular proteins, resulting in a dark blue stained cell.

The Trypan Blue apps on CellDrop Automated Cell Counters enable rapid automated cell counting and viability of cell suspensions stained with Trypan Blue.

### Kit Contents

Kits contain 0.4% Trypan Blue in PBS. The Trypan Blue reagent should be stored at room temperature (15 – 30°C) in an airtight container and does not need to be protected from light.

### Assay Size Trypan Blue Concentration Number of Tests

0.25 mL	0.4%	50
	0.2%	100
1.5 mL	0.4%	300
	0.2%	600

### Sample Volume and Chamber Height

The required sample volume for the CellDrop depends on the height of the measurement chamber, which is set in the counting protocol.

### Standard Magnification (FLi & BF)

#### Gap Height (um) Volume (uL) Minimum Density (cells/mL) Maximum Density (cells/mL)

400	40	7.0E+02	3.1E+06
100	10	2.9E+03	1.3E+07
50	5	5.9E+03	2.5E+07

### Higher Magnification (FLxi & BFx)

#### Gap Height (um) Volume (uL) Minimum Density (cells/mL) Maximum Density (cells/mL)

400	40	4.3E+03	2.6E+07
100	10	1.7E+04	1.0E+08
50	5	3.4E+04	2.1E+08


### Best Practices

- Ensure that the upper and lower chamber surfaces are clean prior to loading sample.
- Lower the arm prior to dispensing sample into the measurement chamber.
- Spin down or filter Trypan Blue dye through a 0.2 µm filter to remove crystallized trypan.
- Mix cells immediately before loading sample and avoid introducing air bubbles.
- Once cells are mixed with trypan blue measure within 5 minutes.
- Follow the image guides to adjust focus and exposure so that unstained cells have bright white centers with a sharp black ring and a sharp transition from light to dark, as shown in Figure 1.
- Allow cells to settle and stop moving across the live preview before pressing the Count button.
- Optimize protocol settings for different cell types. The Default Protocol is a good starting point.

### Sample Prep

Figure 1: Correct focus and exposure settings.

1. Mix cell suspension and Trypan Blue immediately prior to use.



- Optional: Filter Trypan solution through a 0.2  $\mu\text{m}$  filter to remove aggregates and crystals that can form in Trypan solution over time.
- For each sample, mix Trypan and a cell suspension together at the desired ratio and vortex. Refer to the table below for Dilution Factor (DF) guidance examples.

Trypan Volume	Cell Volume	Protocol	Dilution Factor	Recommended Exposure
5 $\mu\text{L}$ 0.4%	5 $\mu\text{L}$	2		Normal
2.5 $\mu\text{L}$ 0.4%	7.5 $\mu\text{L}$	1.33		Low

### Sample Measurement

- With the CellDrop arm in the down position, launch one of the Trypan Blue apps.
- Set sample name, information and protocol as appropriate. If mixing cells and trypan in a ratio other than 1:1, edit the Dilution Factor in the protocol.
- Pipette well-mixed cells + Trypan Blue solution and dispense appropriate sample volume into the measurement chamber, using the groove on the lower sample surface as a pipetting guide.
  - Note:** The volume of sample required depends on the protocol settings for the chamber height. The required volume is displayed on the Count button.
- Adjust exposure and focus according to the image guide.
- Allow cells to settle, then press the Count button.

Refer to [Technical Note 186 – CellDrop Best Practices](#) for additional guidance.

Refer to [denovix.com/sds](https://denovix.com/sds) for safety data sheets for CellDrop Cell Counting Assays.

9-OCT-2024