

# Release LigandTracer Control 2.4

February 23, 2021

## New or Improved

---

- Information of the next phase can be duplicated from MultiDish compartment A/B to compartment C/D.
- Improved handling of files saved to a network folder or external drive: If connection to the folder is lost during a run, LigandTracer Control will try to reconnect each autosave (~ every 10 min). You can manually choose another folder at the end of the run, if the connection has not been reestablished by then.
- Updated handbook

## Fixes

---

- The Experimental Planner is now better suited for runs that uses settings from a previous run (automatic reset of Sent status).

## Known Issues or Limitations

---

- The instrument is disconnected if the computer is put on hibernation or sleep.

